A2657-22 PORTEX NERVE BLOCK -

Smiths Medical ASD, Inc.

APLICARE POVIDONE-IODINE SOLUTION (povidone-iodine solution) solution

[Aplicare, Inc.]

3/4 Ounce Povidone Iodine Packet

Povidone-iodine 10%

Antiseptic

Warnings

Do not use

- if allergic to iodine
- in the eyes

For external use only Ask a doctor before use if injuries are

- deep or puncture wounds
- serious burns

Stop use and ask a doctor if

- redness, irritation, swelling or pain persists or increases
- · infection occurs

Avoid pooling beneath patient

Keep out of reach of children. In case of accidental ingestion, seek professional assistance or consult a poison control center immediately.

Package Label Display Panel



Package Label Display Panel

PORTEX® NERVE BLOCK 22G QUINCKE

PROCEDURAL COMPONENTS:

- 2 Gauze Sponges
- 1 Fenestrated Drape
- 1 Adhesive Bandage
- 1 18G x 1 ½ in. Needle
- 1 25G x 1 1/2 in. Needle
- 1 Needle Stick Pad
- 1 Needle Stop
- 1 Medicine Cup & Cover, 50ml
- 1 Quincke Spinal Needle (22G x 3 ½ in.)
- 1 Plastic Syringe (10ml, Luer Lock)
- 2 Plastic Syringes (5ml, Luer Lock)
- 2 Plastic Syringes (3ml, Luer Lock)

PREP COMPONENTS:

- 1 Povidone-Iodine Solution, ¾ oz.
- 3 Sponge Applicators

WARNINGS:

- As with all plastic wrapped products, do not open or unwr electric discharge may be hazardous, due to risk of explosi
- A needle stick with a contaminated needle may cause infe
- The use of excessive force while placing needles into the s needle to protrude through the bottom of the tray which

PRECAUTIONS:

- Use Aseptic technique.
- To help prevent needle-stick injuries, needles should not be bent. If excessive resistance is met during needle insertion damage may occur. To help avoid needle breakage, do not needle; discard it and complete the procedure with a repla
- After use, place sharps in a suitable sharps container. Dispos in a safe manner according to Centers for Disease Control ar Federal/State/Local regulations (EPA, OSHA) and health ca equivalent.
- Do NOT Resterilize.
- To be used only by individuals familiar with nerve block proce and procedures, refer to standard textbooks.

Smiths design mark and Portex are trademarks of the Smiths The symbol @ indicates the trademark is registered in the U.S. and certain other countries.



Smiths Medical ASD, Inc. Keene, NH 03431 USA Made in USA www.smiths-medical.com





Federal (U.S.A.) law restricts th

the order of a physician. Sterili

