


OXYGEN- oxygen gas
Austin Respiratory Equipment,inc

Disclaimer: This drug has not been found by FDA to be safe and effective, and this labeling has not been approved by FDA. For further information about unapproved drugs, click here.

OXYGEN,
REFRIGERATED LIQUID
USP

UN: 1073



**KEEP OUT OF REACH OF CHILDREN.
DO NOT REMOVE THIS PRODUCT LABEL.**

Produced by air liquefaction.

Directions for use: ↑ Always keep container in upright position ↑. Do not change or force fit connections. Keep oil, grease, and combustibles away. **No smoking in cylinder area.** Use only with equipment cleaned for oxygen service. Store and use with adequate ventilation. Do not get liquid in eyes, on skin, or on clothing. For liquid withdrawal, wear face shield and gloves. Do not drop. Use suitable hand truck for container movement. Avoid spills. Do not walk or roll equipment over spills. Cylinder temperature should not exceed 52°C (125°F). Close valve after each use and when empty. Use a back flow preventive device in the piping. Use in accordance with **Safety data Sheet (SDS).**

FIRST AID: In case of frostbite: Obtain medical treatment immediately.

<u>SIZE</u>	<u>CONTENT</u>
21	16,750 LITERS
31	24,950 LITERS
41	33,500 LITERS

AUSTIN RESPIRATORY EQUIPMENT, INC.
1676 N. Belcher Road
Clearwater, FL 33765
(727) 467-0606

CAS NO:7782-44-7

WARNING

**EXTREMELY COLD, OXIDIZING LIQUID AND GAS UNDER PRESSURE.
VIGOROUSLY ACCELERATES COMBUSTION. COMBUSTIBLES IN
CONTACT WITH LIQUID OXYGEN MAY EXPLODE ON IGNITION OR
IMPACT. CAN CAUSE SEVERE FROSTBITE.**

For emergency use only when administered by properly trained personnel for Oxygen deficiency and resuscitation. For all other medical applications **Rx only** Federal law prohibits dispensing without a prescription. **WARNING: Uninterrupted Use of high concentrations of oxygen over a long duration, without monitoring its Effect on oxygen content of arterial blood, may be harmful. use only with pressure reducing equipment and apparatus designed for oxygen. Do not attempt to use on**

Patients who have stopped breathing, unless used in conjunction with resuscitative equipment. Federal law requires that this Container be Refilled with oxygen U.S.P. only by Establishments Registered as Drug Producer in accordance with the Federal Food, Drug and Cosmetic Act.

KEEP OUT OF REACH OF CHILDREN.

DO NOT REMOVE THIS PRODUCT LABEL.

Produced by Air Liquefaction.

Directions for use: ↑Always keep container in upright position↑. Do not change or force fit connections. Keep oil, grease, and combustibles away. **No smoking in cylinder area.** use only with equipment cleaned for oxygen service. Store and use with adequate ventilation. Do not get liquid in eyes, on skin, or on clothing. For liquid withdrawal, wear face shield and gloves. Do not drop. Use suitable hand truck for container movement.

Avoid spills. Do not walk or roll equipment over spills. Cylinder temperature should not exceed 52°C (125°F).Close valve after each use and when empty. Use a back flow preventive device in the piping. Use in accordance with **Safety data Sheet. (SDS).**

FIRST AID: in case of Frostbite: obtain medical treatment immediately.

<u>SIZE</u>	<u>CONTENT</u>
21	16,750 LITERS
31	24,950 LITERS
41	33,500 LITERS

AUSTIN RESPIRATORY EQUIPMENT, INC.

1676 N. Belcher Road

Clearwater, FL. 33765

(727) 467-0606

OXYGEN				
oxygen gas				
Product Information				
Product Type	HUMAN PRESCRIPTION DRUG	Item Code (Source)	NDC:64540-011(NDC:11390-001)	
Route of Administration	RESPIRATORY (INHALATION)			
Active Ingredient/Active Moiety				
Ingredient Name		Basis of Strength	Strength	
OXYGEN (UNII: S88TT14065) (OXYGEN - UNII:S88TT14065)		OXYGEN	990 mL in 1 L	
Packaging				
#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NDC:64540-011-01	21 L in 1 DEWAR; Type 1: Convenience Kit of Co-Package	07/11/2011	
2	NDC:64540-011-02	41 L in 1 DEWAR; Type 1: Convenience Kit of Co-Package	07/11/2011	
Marketing Information				
Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date	
unapproved medical gas		07/11/2011		

Labeler - Austin Respiratory Equipment,inc (849265483)

Registrant - Austin Respiratory Equipment,inc (849265483)

Establishment

Name	Address	ID/FEI	Business Operations
Austin Respiratory Equipment,inc		849265483	relabel(64540-011)

Revised: 11/2017

Austin Respiratory Equipment,inc