
#### OXYGEN- oxygen gas UnityPoint at Home

Disclaimer: This drug has not been found by FDA to be safe and effective, and this labeling has not been approved by FDA. For further information about unapproved drugs, click here.

-----

#### Oxygen

#### **Principal for Oxygen Product**



Do not handle until all safety precautions have been read and understood. Keep and store away from do thing and other combustible materials. Keep valves and fittings free from grease and oil. Use and store only outdoors or in a well-ventilated place. Wear cold insulating gloves, face shield, and eye protection. In case of fire: Stop leak if safe to do so. Use a back flow preventive device in the piping. Use only with equipment of compatible materials of construction and rated for cylinder pressure. Use only with equipment cleaned for oxygen service. DO NOT change or force fit connections. Avoid spills. Do not walk on or roll equipment over spills. Close valve after each use and when empty. Always keep container in upright position. Read and follow the Safety Data Sheet (SDS) before use.

### FIRST AID:

# IF ON SKIN: <u>Thaw frosted parts with lukewarm water. Do not rub affected</u> area. Get immediate medical advice/attention.

#### DO NOT REMOVE THIS PRODUCT LABEL.

Warning: For emergency use only when administered by properly trained personnel for oxygen deficiency and resuscitation. For all other medical applications, **Rx only.** Uninterrupted use of high concentrations of oxygen over a long duration, without monitoring its effect on oxygen content of arterial blood, may be harmful. Do not attempt to use on patients who have stopped breathing unless used in conjunction with resuscitative equipment. Produced by Air Liquefaction.

Supplied Home Healthcare Equipment 2014

Representative Labeling pursuant to 21 CFR 207.3 (10)
Company name and contact information are printed in this area

OXYGEN,	
REFRIGERATED	
LIQUID USP	
UN 1073	
CONTENTS	LITERS
CAS NO: 7782-44-7	
DANGER:	

REFRIGERATED LIQUEFIED OXIDIZING GAS. ALWAYS
KEEP CONTAINER IN UPRIGHT POSITION. DO NOT
CHANGE OR FORCE FIT CONNECTIONS. MAY CAUSE OR INTENSIFY
FIRE; OXIDIZER, CONTAINS REFRIGERATED GAS: MAY CAUSE
CRYOGENIC BURNS OR INJURY. COMBUSTIBLES IN CONTACT WITH
LIQUID OXYGEN MAY EXPLODE ON IGNITION OR IMPACT.

Do not get in eyes, on skin, or on clothing. Keep reduction valves free from grease and oil. In case of fire, stop leak if safe to do so. Avoid spills. Do not walk or roll equipment over spills. Store and use in a well ventilated place. Use only with equipment of compatible materials of construction and rated for cylinder pressure. Wear cold insulating gloves/face shield/eye protection. Protect from sunlight when ambient temperature exceeds 52 F. Container temperature should not exceed 52 F. Close valve after each use and when empty. Use a back flow preventative device in the piping. Do not handle until all safety precautions have been read and understood. Use only with equipment cleaned for oxygen service. Keep/store away from combustible materials. Do not drop. Use suitable hand truck for container movement.

#### FIRST AID:

IN CASE OF FROSTBITE: Thaw frosted parts with lukewarm water.

Do not rub affected area. Get immediate medical advice/attention.

DO NOT REMOVE THIS PRODUCT LABEL.

For emergency use only when administered by properly trained personnel for oxygen deficiency and resuscitation. For all other medical applications, Rx only. Uninterrupted use of high concentrations of oxygen over a long duration, without monitoring its effect on oxygen content of arterial blood, may be harmful. Use only with pressure reducing equipment and apparatus designed for oxygen. Do not attempt to use on patients who have stopped breathing unless used in conjunction with resuscitative equipment. Produced by Air Liquefaction. ALWAYS KEEP CONTAINER IN UPRIGHT POSITION

## OXYGEN, COMPRESSED USP UN 1072 DANGER: MAY CAUSE OR INTENSIFY FIRE; OXIDIZER. CONTAINS GAS UNDER PRESSURE:

MAY EXPLODE IF HEATED.

Do not handle until all safety precautions have been read and understood Keep and store away from clothing and other combustible materials. Keep valves and fittings free from grease and oil. Use and store only outdoors or in a well-ventilated place. In case of fire: Stop leak if safe to do so. Use a back flow preventive device in the piping. Use only with equipment of compatible materials of construction and rated for cylinder pressure. Use only with equipment cleaned for oxygen service. Open valve slowly. Close valve after each use and when empty. Protect from sunlight when ambient temperature exceeds 52 °C (125 °F). Read and follow the Safety Data Sheet (SDS) before use.

DO NOT REMOVE THIS PRODUCT LABEL.

WARNING: For emergency use only when administered by properly trained personnel for oxygen deficiency and resuscitation. For all other medical applications, Rx only. Uninterrupted use of high concentrations of oxygen over a long duration, without monitoring its effect on oxygen content of arterial blood, may be harmful. Do not attempt to use on patients who have stopped breathing unless used in conjunction with resuscitative equipment. Produced by Air Liquefaction, shooted Home Healthcare Equipment 2014

Contents Size 113 Liters □ M4 165 Liters □ M6 255 Liters □ C 425 Liters □ D 680 Liters DE 7080 Liters □ H Liters ☐ Other

CAS NO: 7782-44-7

Representative Labeling pursuant to 21 CFR 207.3 (10) Company name and contact information are printed in this box

**OXYGEN** 

COMPRESSED

USP

**UN 1072** 

DANGER:

MAY CAUSE OR INTENSIFY FIRE; OXIDIZER. CONTAINS GAS UNDER PRESSURE; MAY EXPLODE IF HEATED.

Do not handle until all safety precautions have been read and understood. Keep and store away from clothing and other combustible materials. Keep valves and fittings free from grease and oil. Use and store only outdoors or in a well-ventilated place. In case of fire: Stop leak if safe to do so. Use a back flow preventive device in the piping. Use only with equipment of compatible materials of construction and rated for cylinder pressure. Use only with equipment cleaned for oxygen service. Open valve slowly. Close valve after each use and when empty. Protect from sunlight when ambient temperature exceeds 52 °c (125 °F). Read and follow the Safety Data Sheet (SDS) before use.

CAS NO: 7782-44-7

DO NOT REMOVE THIS PRODUCT LABEL.

WARNING: For emergency use only when administered by properly trained personnel for oxygen deficiency

and resuscitation. For all other medical applications, Rx only.

Uninterrupted use of high concentrations of oxygen over a long duration, without monitoring its effect on oxygen

content of arterial blood, may be harmful. Do not attempt to use on patients who have stopped breathing unless used in conjunction with resuscitative equipment. Produced by Air Liquefaction.



Representative Labeling pursuant to 21 CFR 207.3 (10)

#### **OXYGEN**

oxygen gas

Product Information			
Product Type	HUMAN PRESCRIPTION DRUG	Item Code (Source)	NDC:66698-001
Route of Administration	RESPIRATORY (INHALATION)		

Active Ingredient/Active Moiety			
Ingredient Name	<b>Basis of Strength</b>	Strength	
<b>OXYGEN</b> (UNII: S88TT14065) (OXYGEN - UNII:S88TT14065)	OXYGEN	99 L in 100 L	

Packaging				
#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NDC:66698- 001-02	164 L in 1 CYLINDER; Type 0: Not a Combination Product	10/01/2014	
7	NDC:66698-	165 L in 1 CYLINDER; Type 0: Not a Combination	10/01/2014	

4	001-01	Product	10/01/2014	
3	NDC:66698- 001-03	255 L in 1 CYLINDER; Type 0: Not a Combination Product	10/01/2014	
4	NDC:66698- 001-04	425 L in 1 CYLINDER; Type 0: Not a Combination Product	10/01/2014	
5	NDC:66698- 001-05	680 L in 1 CYLINDER; Type 0: Not a Combination Product	10/01/2014	
6	NDC:66698- 001-06	1738 L in 1 CYLINDER; Type 0: Not a Combination Product	10/01/2014	
7	NDC:66698- 001-07	3455 L in 1 CYLINDER; Type 0: Not a Combination Product	10/01/2014	
8	NDC:66698- 001-08	7080 L in 1 CYLINDER; Type 0: Not a Combination Product	10/01/2014	
9	NDC:66698- 001-09	20 L in 1 DEWAR; Type 0: Not a Combination Product	10/01/2014	
10	NDC:66698- 001-10	40 L in 1 DEWAR; Type 0: Not a Combination Product	10/01/2014	
11	NDC:66698- 001-11	60 L in 1 DEWAR; Type 0: Not a Combination Product	10/01/2014	
12	NDC:66698- 001-12	180 L in 1 DEWAR; Type 0: Not a Combination Product	10/01/2014	
13	NDC:66698- 001-13	719 L in 1 DEWAR; Type 0: Not a Combination Product	10/01/2014	
14	NDC:66698- 001-14	21 L in 1 DEWAR; Type 0: Not a Combination Product	10/01/2014	
15	NDC:66698- 001-15	41 L in 1 DEWAR; Type 0: Not a Combination Product	10/01/2014	
16	NDC:66698- 001-16	46 L in 1 DEWAR; Type 0: Not a Combination Product	10/01/2014	

Marketing Information			
Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
unapproved medical gas		10/01/2014	

## Labeler - UnityPoint at Home (151070955)

Establishment				
Name	Address	ID/FEI	<b>Business Operations</b>	
UnityPoint at Home		061943101	manufacture(66698-001)	

Revised: 10/2023 UnityPoint at Home