

**TD GOLD GLUCOSE CONTROL SOLUTION LEVEL 2 SET- glucose dehydrogenase, glucose**  
**EPS Bio Technology Corp.**

-----  
**Home Use Medical Device- TD Gold Glucose Control Solution Level 2**

**TD Gold Glucose Control Solution Level 2 Box**



**TD Gold Glucose Control Solution Level 2 label**

# SIZE: 55 X 30 mm

# TD Gold

## Glucose Control Solutions

# Level 2 4ml

Recap the bottle immediately after use.  
Discard test strips in 90 days after first opening date.

Open on /

LOT

EXP

- TD Gold Glucose Control Solution is for in vitro diagnostic use only.
- Rock the control solution bottle back and forth several times before use.
- DO NOT freeze.
- Use the control solution only for testing with TD Gold Meters and Test Strips.
- Store the TD Gold Glucose Control Solution between 36-86°F (2-30°C).

Manufactured for Total Diabetes, LLC  
Telephone Assistance / 855-525-8332

Manufacturer / EPS Bio Technology Corp.  
PN 706400019\_1249\_02

## TD GOLD GLUCOSE CONTROL SOLUTION LEVEL 2 SET

glucose dehydrogenase, glucose

### Product Information

Product Type	MEDICAL DEVICE	Item Code (Source)	NHRIC:47884-184
--------------	----------------	--------------------	-----------------

### Product Characteristics

number of times usable	80
------------------------	----

### Packaging

#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NHRIC:47884-184-00	1 in 1 PACKAGE		

### Marketing Information

Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
premarket notification	K101037	10/01/2013	

**Labeler** - EPS Bio Technology Corp. (658620117)

