

EYESALINE CONCENTRATE- water liquid
Honeywell Safety Products USA, Inc

0498-0601: Eyesaline Concentrate

Active Ingredient

Purified water 93%

Purpose

Eyewash

Uses

- For flushing or irrigating the eye to remove loose foreign material, air pollutants, or chlorinated water.

Warnings

For external use only- Obtain immediate medical treatment for all open wounds in or near the eyes.

Do not use

- in concentrated form
- while wearing contact lenses
- if solution changes color or becomes cloudy
- for injection
- in intraocular surgery
- internally
- if eye is lacerated or object is embedded
- if seal is broken

Stop use and consult a doctor if you have

- changes in vision
- continued redness, irritation or pain, or if the condition worsens or persists

Keep out of reach of children

If swallowed, get medical help or contact a Poison Control Center right away.

Directions

- Concentrate must be diluted (1) part concentrate to (10) parts potable water before use
- Prepare sufficient quantity of constituted eyewash solution in advance of an emergency

- Mix concentrate with potable water to prepare eyewash solution as follows:
- Determine quantity of eyewash solution required
- Follow sanitary procedures
- Wear protective eyewear and gloves
- Pour 1/3 of required amount of potable water into mixing container or eyewash unit
- Add 1/3 of required amount of concentrate into water
- Agitate or mix thoroughly
- Repeat until required amount of solution is prepared
- Properly dispose of empty bottle after use

Do not reuse.

- Replace solution no longer than (6) months after mixing

Other information

- Do not freeze

Inactive ingredients

benzalkonium chloride, edetate disodium, sodium chloride, sodium phosphate dibasic, sodium phosphate monobasic

Questions?

Call **1-800-430-5490**

Sperian Eye and Face Protection Inc

(a Honeywell Company)

825 East Highway 151

Platteville, WI 53818 USA

Package Label

Honeywell

DO NOT USE IF SEAL IS BROKEN OR TAMPERED WITH
NO USAR SI EL SELLO ÉSTA ROTO O HA SIDO ALTERADO
NE PAS UTILISER SI LE SCEAU EST BRISE OU A ÊTÈ ALTERÉ

eyesaline®

CONCENTRATE CONCENTRADO CONCENTRÉ

When mixed according to directions, Eyesaline Concentrate and water will provide a buffered, preserved saline solution for flushing or irrigating the eyes.

Cuando se mezcla según las instrucciones Concentrado Eyesaline y agua, se obtiene una solución salina tamponada y preservada para enjuagar o irrigar los ojos.

Lorsque le concentré Eyesaline et l'eau sont mélangés selon les instructions, ce mélange fournira une solution saline tamponnée, préservée pour rincer ou irriguer les yeux.

70 fl. oz. (2.07 L)

REORDER / NUEVO PEDIDO / RÉAPPROVISIONNEMENT #32-000509-0000

#32-004622 Rev. E

EYESALINE CONCENTRATE

water liquid

Product Information

Product Type	HUMAN OTC DRUG	Item Code (Source)	NDC:0498-0601
Route of Administration	OPHTHALMIC		

Active Ingredient/Active Moiety

Ingredient Name	Basis of Strength	Strength
WATER (UNII: 059QF0KO0R) (WATER - UNII:059QF0KO0R)	WATER	92.2 L in 100 L

Inactive Ingredients

Ingredient Name	Strength
SODIUM PHOSPHATE, MONOBASIC, MONOHYDRATE (UNII: 593YOG76RN)	
SODIUM PHOSPHATE, DIBASIC, ANHYDROUS (UNII: 22ADO53M6F)	
SODIUM CHLORIDE (UNII: 451W47IQ8X)	
EDETATE DISODIUM (UNII: 7FLD91C86K)	
BENZALKONIUM CHLORIDE (UNII: F5UM2KM3W7)	

Packaging

#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NDC:0498-0601-18	2.07 L in 1 BOTTLE, PLASTIC; Type 0: Not a Combination Product	09/01/2020	
2	NDC:0498-0601-19	5.32 L in 1 BOTTLE, PLASTIC; Type 0: Not a Combination Product	09/01/2020	

Marketing Information

Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
OTC Monograph Drug	M018	09/01/2020	

Labeler - Honeywell Safety Products USA, Inc (118768815)

Establishment

Name	Address	ID/FEI	Business Operations
Keystone Industries		014769301	manufacture(0498-0601)

Revised: 1/2024

Honeywell Safety Products USA, Inc