

**ACD-A- anticoagulant citrate dextrose solution a solution**  
**Haemonetics Manufacturing Inc**

-----  
**Anticoagulant Citrate Dextrose Solution, USP (ACD), Solution A**

**DESCRIPTION**

Re-Order Product Code: 796-08

750 mL volume. Single use container.

Sterile, nonpyrogenic fluid path. Sterilized by steam.

Store at room temperature.

**INDICATIONS AND USAGE**

For use with automated pheresis equipment only. Use according to equipment manufacturer's instructions. This unit should be used for pheresis only.

**CONTRAINDICATIONS**

Not for direct intravenous infusion.

**WARNINGS**

Avoid excessive heat. Protect from freezing.

**PRECAUTIONS**

Do not use unless solution is clear and no leaks detected.

**General**

Rx only. Discard unused portion.

Recommended Storage: Store at room temperature (25°C / 77°F).

**HOW SUPPLIED**

750 mL volume bag.

Each 100 ml contains:

0.73 g citric acid (anhydrous), USP

2.20 g sodium citrate (dihydrate), USP

2.45 g dextrose (monohydrate), USP

**PACKAGE/LABEL PRINCIPAL DISPLAY PANEL**

Hot Stamped Bag

# HAEMONETICS®

Anticoagulant Citrate Dextrose  
Solution, USP (ACD)  
Solution A—750 ml

- |    |  |    |
|----|--|----|
| 1— | Each 100 ml contains:<br>0.73 g citric acid (anhydrous), USP<br>2.20 g sodium citrate (dihydrate), USP<br>2.45 g dextrose (monohydrate), USP   | —1 |
| 2— | Sterile, nonpyrogenic fluid path. Sterilized<br>by steam. <b>Rx Only.</b>  | —2 |
| 3— | Single use container. Discard unused<br>portion.   | —3 |
| 4— | Not for direct intravenous infusion.<br>Do not use unless solution is clear and<br>no leaks are detected.  | —4 |
| 5— | Recommended storage: Room temperature<br>(25 °C / 77 °F). Avoid excessive heat.<br>Protect from freezing.  | —5 |
|    | Code 796-08                      DIN 02208105  |    |
| 6— | Manufactured for: Haemonetics Corporation<br>400 Wood Road<br>Braintree, MA 02184, USA<br>By: Haemonetics Manufacturing Inc.<br>1650 Industrial Park Road<br>Covington, LA 70122, USA<br><br>143796086 | —6 |
| 7— |  | —7 |

## Envelope Label

### HAEMONETICS®

LOT #####

12 Units                      YYYYY-MM #

Anticoagulant Citrate Dextrose Solution,  
USP (ACD), Solution A - 750 ml

**Rx only.** Sterile, nonpyrogenic fluid path. Sterilized by  
steam. Store at room temperature. Avoid excessive heat.  
Protect from freezing.

Manufactured for: Haemonetics Corporation  
400 Wood Road, Braintree, MA 02184, USA  
By: Haemonetics Manufacturing Inc.  
1650 Industrial Park Road  
Covington, LA 70122, USA

152790284

Code 796-08

[                      Insert Barcode Here                      ]

02110088/0210005/01001701/YYMMDD10#####

## References

HAEMONETICS and THE Blood Management Company are trademarks of Haemonetics Corporation.

Manufactured for:

Haemonetics Corporation

400 Wood Road

Braintree, MA 02184, USA

By: Haemonetics Manufacturing Inc.

1630 Industrial Park Street  
Covina, CA 91722, USA

## ACD-A

anticoagulant citrate dextrose solution a solution

### Product Information

|                                |                         |                           |               |
|--------------------------------|-------------------------|---------------------------|---------------|
| <b>Product Type</b>            | HUMAN PRESCRIPTION DRUG | <b>Item Code (Source)</b> | NDC:53157-796 |
| <b>Route of Administration</b> | INTRAVENOUS             |                           |               |

### Active Ingredient/Active Moiety

| Ingredient Name  | Basis of Strength     | Strength            |
|--|-----------------------|---------------------|
| DEXTROSE MONOHYDRATE (UNII: LX22YL083G) (ANHYDROUS DEXTROSE - UNII:5SL0G7R0OK)       | DEXTROSE MONOHYDRATE  | 2.45 g<br>in 100 mL |
| SODIUM CITRATE (UNII: 1Q73Q2JULR) (ANHYDROUS CITRIC ACID - UNII:XF417D3PSL)          | SODIUM CITRATE        | 2.20 g<br>in 100 mL |
| CITRIC ACID MONOHYDRATE (UNII: 2968PHW8QP) (ANHYDROUS CITRIC ACID - UNII:XF417D3PSL) | ANHYDROUS CITRIC ACID | 0.73 g<br>in 100 mL |

### Packaging

| # | Item Code        | Package Description                                | Marketing Start Date | Marketing End Date |
|---|------------------|--|----------------------|--------------------|
| 1 | NDC:53157-796-08 | 12 in 1 CARTON                                     |                      |                    |
| 1 |                  | 750 mL in 1 BAG; Type 0: Not a Combination Product |                      |                    |

### Marketing Information

| Marketing Category | Application Number or Monograph Citation | Marketing Start Date | Marketing End Date |
|--------------------|--|----------------------|--------------------|
| NDA                | NDA710497                                | 11/06/1987           |                    |

**Labeler** - Haemonetics Manufacturing Inc (078598396)

**Registrant** - Haemonetics Manufacturing Inc (078598396)

### Establishment

| Name                          | Address | ID/FEI    | Business Operations    |
|-------------------------------|---------|-----------|------------------------|
| Haemonetics Manufacturing Inc |         | 078598396 | MANUFACTURE(53157-796) |