

**HAEMONETICS ANTICOAGULANT SODIUM CITRATE- trisodium citrate dihydrate solution**

**Haemonetics Corporation**

-----

**Haemonetics Anticoagulant Sodium Citrate**

**HAEMONETICS ANTICOAGULANT SODIUM CITRATE 4% w/v SOLUTION USP**

**Rx only                      250 mL**

**Intended for use only with automated apheresis devices.**

**Each 100 mL contains:** Sodium Citrate (Dihydrate), USP 4.0g  
(pH adjusted with Citric Acid, Monohydrate, USP)

**CAUTION:**

Not for direct intravenous infusion. The pouch is a moisture barrier. Do not use unless solution is clear and no leaks detected. Single use container. Discard unused portion.

STERILE, nonpyrogenic fluid path.

**RECOMMENDED STORAGE:**

Room temperature (25°C/77°F). Avoid excess heat. Protect from freezing.

**Product code**

Product Code 420A  
117881-00, Rev. AA-XXX

Haemonetics Corporation  
400 Wood Road  
Braintree, MA 02184 USA

**Product Labeling**

**HAEMONETICS<sup>®</sup>**  
**ANTICOAGULANT**  
**SODIUM CITRATE 4% w/v**  
**SOLUTION U.S.P.**

**250 mL**

**Intended for use only with automated apheresis devices**

**Each 100 mL contains:**

Sodium Citrate (Dihydrate), USP      4.0 g  
(pH adjusted with Citric Acid, Monohydrate, USP)

**CAUTION:** Not for direct intravenous infusion. The pouch is a moisture barrier. Do not use unless solution is clear and no leaks detected. Single use container. Discard unused portion.

**STERILE**, nonpyrogenic fluid path.

Rx only.

**RECOMMENDED STORAGE:** Room temperature (25°C/77°F).  
Avoid excessive heat. Protect from freezing.

Product code 420A

112534-00, Rev. AA / M040, Rev. 01/12

Haemonetics Corporation  
400 Wood Road  
Braintree, MA 02184 USA  
Made in Mexico

Lot No.

Exp. Date

A rectangular box with a thin black border, intended for handwritten entry of the Lot No. and Exp. Date. The box is currently empty.



# HAEMONETICS®

## 420A (Romania) Product Bag



## 420A (Romania) Product Barcoding



**HAEMONETICS**®

**ANTICOAGULANT  
SODIUM CITRATE  
4% w/v SOLUTION USP**

**Rx only                      250 mL**

**Intended for use only with automated apheresis devices.**

**Each 100 mL contains:**

**Sodium Citrate (Dihydrate), USP 4.0g**

**(pH adjusted with Citric Acid, Monohydrate, USP)**

**CAUTION: Not for direct intravenous infusion.**

**The pouch is a moisture barrier. Do not use unless solution is clear and no leaks detected. Single use container. Discard unused portion.**

**STERILE, nonpyrogenic fluid path.**

**RECOMMENDED STORAGE: Room temperature (25°C/77°F).**

**Avoid excess heat. Protect from freezing.**

**Product code 420A  
117881-00, Rev. AA-XXX**

**Haemonetics Corporation  
400 Wood Road  
Braintree, MA 02184, USA**

<b>HAEMONETICS ANTICOAGULANT SODIUM CITRATE</b>			
trisodium citrate dihydrate solution			
<b>Product Information</b>			
<b>Product Type</b>	HUMAN PRESCRIPTION DRUG	<b>Item Code (Source)</b>	NDC:57826-420
<b>Route of Administration</b>	EXTRACORPOREAL		
<b>Active Ingredient/Active Moiety</b>			
	<b>Ingredient Name</b>	<b>Basis of Strength</b>	<b>Strength</b>
	TRISODIUM CITRATE DIHYDRATE (UNII: B22547B95K) (ANHYDROUS CITRIC ACID - UNII:XF417D3PSL)	ANHYDROUS CITRIC ACID	40 mg in 1 mL
<b>Inactive Ingredients</b>			

Ingredient Name	Strength
WATER (UNII: 059QF0KO0R)	

**Packaging**

#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NDC:57826-420-02	250 mL in 1 BAG; Type 0: Not a Combination Product		

**Marketing Information**

Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
NDA	BN980123	01/10/2013	

**Labeler** - Haemonetics Corporation (057827420)