

**SINBIO SANITIZER S-300- didecyldimonium chloride and benzalkonium chloride liquid liquid
SINBIO, S.A. DE C.V.**

Disclaimer: Most OTC drugs are not reviewed and approved by FDA, however they may be marketed if they comply with applicable regulations and policies. FDA has not evaluated whether this product complies.

DOSAGE & ADMINISTRATION SECTION

SANITIZER S-300

VIRUCIDAL + GERMICIDAL

Hand and Surface Sanitizer

Hospitals | Restaurants | Offices | Home



Active Ingredients:

Alkyl Dymethyl Benzyl
Ammomium Chloride C12 - C16 (32%)
Benzalkonium Chloride (10%)

Purpose

Antiseptic

Indications & usage:

Use it directly from the aplicator on the hands, body, surfaces or objects such as cellphones, keyboards, door handles, any utencil or tool. Works as a better alternative than an alcohol based hand sanitizer. Drying time is 20 seconds at 25°C. Recommended for repeated use.

Other Information:

Form: Liquid
Color: Cristal Blue
pH: Neutral
Solubility: 100%
Temperature: 15-40°C
Specific Weight: 1.084 kg/L
Scent: Citrus

Inactive Ingredients:

Water, Octy decyl dimethyl ammonium chloride, Dioctyl dimethyl ammonium chloride, Didecyl dimethyl ammonium chloride, D-Glucopyranose olimers glycosides, Propan 1, 2, 3 - trial, Triarylmethane.

Warnings:

Keep out of children reach. If swallowed, get medical help or contact a Poison Control Center right away. When using, keep out of eyes. In case of eye contact, rinse eyes thoughly with water. EXTERNAL USE ONLY. DO NOT INGEST THIS PRODUCT. May cause eye irritation.



Dosage & Administration:

Apply sanitizer to palm. Rub the product over the surfaces of both hands until dry. To create foam, product must be dispensed from a foaming pump. Product can also be dispensed in a liquid consistency using a spray bottle.



Eliminates 99.99%
Virus & Bacteria

sinbio  /sinbiomid

(999) 909 4282
(999) 981 4627

MÉRIDA, YUCATÁN, MX
www.sinbio.com.mx
contacto@sinbio.com.mx

Manufactured by:
Sinbio, S. A. de C. V.
FDA 79370

NET CONT. :

1.32gal/5L

LOT AND EXP DATE:

79370-001-01

SANITIZER S-300

VIRUCIDAL + GERMICIDAL

Hand and Surface Sanitizer

Hospitals | Restaurants | Offices | Home



Active Ingredients:

Alkyl Dymethyl Benzyl
Ammomium Chloride C12 - C16 (32%)
Benzalkonium Chloride (10%)

Purpose

Antiseptic

Indications & usage:

Use it directly from the aplicator on the hands, body, surfaces or objects such as cellphones, keyboards, door handles, any utencil or tool. Works as a better alternative than an alcohol based hand sanitizer. Drying time is 20 seconds at 25°C. Recommended for repeated use.

Other Information:

Form: Liquid
Color: Cristal Blue
pH: Neutral
Solubility: 100%
Temperature: 15-40°C
Specific Weight: 1.084 kg/L
Scent: Citrus

Inactive Ingredients:

Water, Octy decyl dimethyl ammonium chloride, Dioctyl dimethyl ammonium chloride, Didecyl dimethyl ammonium chloride, D-Glucopyranose olimers glycosides, Propan 1, 2, 3 - trial, Triarylmethane.

Warnings:

Keep out of children reach. If swallowed, get medical help or contact a Poison Control Center right away. When using, keep out of eyes. In case of eye contact, rinse eyes thoughly with water. EXTERNAL USE ONLY. DO NOT INGEST THIS PRODUCT. May cause eye irritation.



Dosage & Administration:

Apply sanitizer to palm. Rub the product over the surfaces of both hands until dry. To create foam, product must be dispensed from a foaming pump. Product can also be dispensed in a liquid consistency using a spray bottle.



Eliminates 99.99%
Virus & Bacteria

sinbio  /sinbiomid

(999) 909 4282
(999) 981 4627

MÉRIDA, YUCATÁN, MX
www.sinbio.com.mx
contacto@sinbio.com.mx

Manufactured by:
Sinbio, S. A. de C. V.
FDA 79370

NET CONT. :

1.32gal/5L

LOT AND EXP DATE:

79370-001-01

SANITIZER S-300

VIRUCIDAL + GERMICIDAL

Hand and Surface Sanitizer

Hospitals | Restaurants | Offices | Home



Active Ingredients:

Alkyl Dymethyl Benzyl
Ammomium Chloride C12 - C16 (32%)
Benzalkonium Chloride (10%)

Purpose

Antiseptic

Indications & usage:

Use it directly from the aplicator on the hands, body, surfaces or objects such as cellphones, keyboards, door handles, any utencil or tool. Works as a better alternative than an alcohol based hand sanitizer. Drying time is 20 seconds at 25°C. Recommended for repeated use.

Other Information:

Form: Liquid
Color: Cristal Blue
pH: Neutral
Solubility: 100%
Temperature: 15-40°C
Specific Weight: 1.084 kg/L
Scent: Citrus

Inactive Ingredients:

Water, Octy decyl dimethyl ammonium chloride, Dioctyl dimethyl ammonium chloride, Didecyl dimethyl ammonium chloride, D-Glucopyranose olimers glycosides, Propan 1, 2, 3 - trial, Triarylmethane.

Warnings:

Keep out of children reach. If swallowed, get medical help or contact a Poison Control Center right away. When using, keep out of eyes. In case of eye contact, rinse eyes thoughly with water. EXTERNAL USE ONLY. DO NOT INGEST THIS PRODUCT. May cause eye irritation.



Dosage & Administration:

Apply sanitizer to palm. Rub the product over the surfaces of both hands until dry. To create foam, product must be dispensed from a foaming pump. Product can also be dispensed in a liquid consistency using a spray bottle.



Eliminates 99.99%
Virus & Bacteria

sinbio  /sinbiomid

(999) 909 4282
(999) 981 4627

MÉRIDA, YUCATÁN, MX
www.sinbio.com.mx
contacto@sinbio.com.mx

Manufactured by:
Sinbio, S. A. de C. V.
FDA 79370

NET CONT. :

1.32gal/5L

LOT AND EXP DATE:

79370-001-01

SANITIZER S-300

VIRUCIDAL + GERMICIDAL

Hand and Surface Sanitizer

Hospitals | Restaurants | Offices | Home



Active Ingredients:

Alkyl Dymethyl Benzyl
Ammomium Chloride C12 - C16 (32%)
Benzalkonium Chloride (10%)

Purpose

Antiseptic

Indications & usage:

Use it directly from the aplicator on the hands, body, surfaces or objects such as cellphones, keyboards, door handles, any utencil or tool. Works as a better alternative than an alcohol based hand sanitizer. Drying time is 20 seconds at 25°C. Recommended for repeated use.

Other Information:

Form: Liquid
Color: Cristal Blue
pH: Neutral
Solubility: 100%
Temperature: 15-40°C
Specific Weight: 1.084 kg/L
Scent: Citrus

Inactive Ingredients:

Water, Octy decyl dimethyl ammonium chloride, Dioctyl dimethyl ammonium chloride, Didecyl dimethyl ammonium chloride, D-Glucopyranose olimers glycosides, Propan 1, 2, 3 - trial, Triarylmethane.

Warnings:

Keep out of children reach. If swallowed, get medical help or contact a Poison Control Center right away. When using, keep out of eyes. In case of eye contact, rinse eyes thoughly with water. EXTERNAL USE ONLY. DO NOT INGEST THIS PRODUCT. May cause eye irritation.



Dosage & Administration:

Apply sanitizer to palm. Rub the product over the surfaces of both hands until dry. To create foam, product must be dispensed from a foaming pump. Product can also be dispensed in a liquid consistency using a spray bottle.



Eliminates 99.99%
Virus & Bacteria

sinbio  /sinbiomid

(999) 909 4282
(999) 981 4627

MÉRIDA, YUCATÁN, MX
www.sinbio.com.mx
contacto@sinbio.com.mx

Manufactured by:
Sinbio, S. A. de C. V.
FDA 79370

NET CONT. :

1.32gal/5L

LOT AND EXP DATE:

79370-001-01

OTC - KEEP OUT OF REACH OF CHILDREN SECTION

SANITIZER S-300

VIRUCIDAL + GERMICIDAL

Hand and Surface Sanitizer

Hospitals | Restaurants | Offices | Home



Active Ingredients:

Alkyl Dymethyl Benzyl
Ammomium Chloride C12 - C16 (32%)
Benzalkonium Chloride (10%)

Purpose

Antiseptic

Indications & usage:

Use it directly from the aplicator on the hands, body, surfaces or objects such as cellphones, keyboards, door handles, any utencil or tool. Works as a better alternative than an alcohol based hand sanitizer. Drying time is 20 seconds at 25°C. Recommended for repeated use.

Other Information:

Form: Liquid
Color: Cristal Blue
pH: Neutral
Solubility: 100%
Temperature: 15-40°C
Specific Weight: 1.084 kg/L
Scent: Citrus

Inactive Ingredients:

Water, Octy decyl dimethyl ammonium chloride, Dioctyl dimethyl ammonium chloride, Didecyl dimethyl ammonium chloride, D-Glucopyranose olimers glycosides, Propan 1, 2, 3 - trial, Triarylmethane.

Warnings:

Keep out of children reach. If swallowed, get medical help or contact a Poison Control Center right away. When using, keep out of eyes. In case of eye contact, rinse eyes thoughly with water. EXTERNAL USE ONLY. DO NOT INGEST THIS PRODUCT. May cause eye irritation.



Dosage & Administration:

Apply sanitizer to palm. Rub the product over the surfaces of both hands until dry. To create foam, product must be dispensed from a foaming pump. Product can also be dispensed in a liquid consistency using a spray bottle.



Eliminates 99.99%
Virus & Bacteria

sinbio  /sinbiomid

(999) 909 4282
(999) 981 4627

MÉRIDA, YUCATÁN, MX
www.sinbio.com.mx
contacto@sinbio.com.mx

Manufactured by:
Sinbio, S. A. de C. V.
FDA 79370

NET CONT. :

1.32gal/5L

LOT AND EXP DATE:

79370-001-01

SANITIZER S-300

VIRUCIDAL + GERMICIDAL

Hand and Surface Sanitizer

Hospitals | Restaurants | Offices | Home



Active Ingredients:

Alkyl Dymethyl Benzyl
Ammomium Chloride C12 - C16 (32%)
Benzalkonium Chloride (10%)

Purpose

Antiseptic

Indications & usage:

Use it directly from the aplicator on the hands, body, surfaces or objects such as cellphones, keyboards, door handles, any utencil or tool. Works as a better alternative than an alcohol based hand sanitizer. Drying time is 20 seconds at 25°C. Recommended for repeated use.

Other Information:

Form: Liquid
Color: Cristal Blue
pH: Neutral
Solubility: 100%
Temperature: 15-40°C
Specific Weight: 1.084 kg/L
Scent: Citrus

Inactive Ingredients:

Water, Octy decyl dimethyl ammonium chloride, Dioctyl dimethyl ammonium chloride, Didecyl dimethyl ammonium chloride, D-Glucopyranose olimers glycosides, Propan 1, 2, 3 - trial, Triarylmethane.

Warnings:

Keep out of children reach. If swallowed, get medical help or contact a Poison Control Center right away. When using, keep out of eyes. In case of eye contact, rinse eyes thoughly with water. EXTERNAL USE ONLY. DO NOT INGEST THIS PRODUCT. May cause eye irritation.



Dosage & Administration:

Apply sanitizer to palm. Rub the product over the surfaces of both hands until dry. To create foam, product must be dispensed from a foaming pump. Product can also be dispensed in a liquid consistency using a spray bottle.



Eliminates 99.99%
Virus & Bacteria

sinbio /sinbiomid

(999) 909 4282
(999) 981 4627

MÉRIDA, YUCATÁN, MX
www.sinbio.com.mx
contacto@sinbio.com.mx

Manufactured by:
Sinbio, S. A. de C. V.
FDA 79370

NET CONT. :

1.32gal/5L

LOT AND EXP DATE:

79370-001-01

SANITIZER S-300

VIRUCIDAL + GERMICIDAL

Hand and Surface Sanitizer

Hospitals | Restaurants | Offices | Home



Active Ingredients:

Alkyl Dymethyl Benzyl
Ammomium Chloride C12 - C16 (32%)
Benzalkonium Chloride (10%)

Purpose

Antiseptic

Indications & usage:

Use it directly from the aplicator on the hands, body, surfaces or objects such as cellphones, keyboards, door handles, any utencil or tool. Works as a better alternative than an alcohol based hand sanitizer. Drying time is 20 seconds at 25°C. Recommended for repeated use.

Other Information:

Form: Liquid
Color: Cristal Blue
pH: Neutral
Solubility: 100%
Temperature: 15-40°C
Specific Weight: 1.084 kg/L
Scent: Citrus

Inactive Ingredients:

Water, Octy decyl dimethyl ammonium chloride, Dioctyl dimethyl ammonium chloride, Didecyl dimethyl ammonium chloride, D-Glucopyranose olimers glycosides, Propan 1, 2, 3 - trial, Triarylmethane.

Warnings:

Keep out of children reach. If swallowed, get medical help or contact a Poison Control Center right away. When using, keep out of eyes. In case of eye contact, rinse eyes thoughly with water. EXTERNAL USE ONLY. DO NOT INGEST THIS PRODUCT. May cause eye irritation.



Dosage & Administration:

Apply sanitizer to palm. Rub the product over the surfaces of both hands until dry. To create foam, product must be dispensed from a foaming pump. Product can also be dispensed in a liquid consistency using a spray bottle.



Eliminates 99.99%
Virus & Bacteria

sinbio  /sinbiomid

(999) 909 4282
(999) 981 4627

MÉRIDA, YUCATÁN, MX
www.sinbio.com.mx
contacto@sinbio.com.mx

Manufactured by:
Sinbio, S. A. de C. V.
FDA 79370

NET CONT. :

1.32gal/5L

LOT AND EXP DATE:

79370-001-01

SANITIZER S-300

VIRUCIDAL + GERMICIDAL

Hand and Surface Sanitizer

Hospitals | Restaurants | Offices | Home



Active Ingredients:

Alkyl Dymethyl Benzyl
Ammomium Chloride C12 - C16 (32%)
Benzalkonium Chloride (10%)

Purpose

Antiseptic

Indications & usage:

Use it directly from the aplicator on the hands, body, surfaces or objects such as cellphones, keyboards, door handles, any utencil or tool. Works as a better alternative than an alcohol based hand sanitizer. Drying time is 20 seconds at 25°C. Recommended for repeated use.

Other Information:

Form: Liquid
Color: Cristal Blue
pH: Neutral
Solubility: 100%
Temperature: 15-40°C
Specific Weight: 1.084 kg/L
Scent: Citrus

Inactive Ingredients:

Water, Octy decyl dimethyl ammonium chloride, Dioctyl dimethyl ammonium chloride, Didecyl dimethyl ammonium chloride, D-Glucopyranose olimers glycosides, Propan 1, 2, 3 - trial, Triarylmethane.

Warnings:

Keep out of children reach. If swallowed, get medical help or contact a Poison Control Center right away. When using, keep out of eyes. In case of eye contact, rinse eyes thoughly with water. EXTERNAL USE ONLY. DO NOT INGEST THIS PRODUCT. May cause eye irritation.



Dosage & Administration:

Apply sanitizer to palm. Rub the product over the surfaces of both hands until dry. To create foam, product must be dispensed from a foaming pump. Product can also be dispensed in a liquid consistency using a spray bottle.



Eliminates 99.99%
Virus & Bacteria

sinbio  /sinbiomid

(999) 909 4282
(999) 981 4627

MÉRIDA, YUCATÁN, MX
www.sinbio.com.mx
contacto@sinbio.com.mx

Manufactured by:
Sinbio, S. A. de C. V.
FDA 79370

NET CONT. :

1.32gal/5L

LOT AND EXP DATE:

79370-001-01

SINBIO SANITIZER S-300

didecylmonium chloride and benzalkonium chloride liquid liquid

Product Information

Product Type	HUMAN OTC DRUG	Item Code (Source)	NDC:79370-001
Route of Administration	TOPICAL		

Active Ingredient/Active Moiety

Ingredient Name	Basis of Strength	Strength
N-ALKYL DIMETHYL BENZYL AMMONIUM CHLORIDE (C12-C18) (UNII: 9U1Q4T4ZYS) (DIDECYLDIMONIUM - UNII:Z7F472XQPA)	N-ALKYL DIMETHYL BENZYL AMMONIUM CHLORIDE (C12-C18)	0.32 mg in 0.001 L
BENZALKONIUM CHLORIDE (UNII: F5UM2KM3W7) (BENZALKONIUM - UNII:7N6JUD5X6Y)	BENZALKONIUM CHLORIDE	0.1 mg in 0.001 L

Inactive Ingredients

Ingredient Name	Strength
OCTYLDIMETHYLAMMONIUM CHLORIDE (UNII: VBX320461J)	
BENZYLDIMETHYLOCTYLAMMONIUM CHLORIDE (UNII: OC19CZY82T)	
WATER (UNII: 059QF0KO0R)	
DIDECYLDIMONIUM CHLORIDE (UNII: JXN40O9Y9B)	
ALCOHOL (UNII: 3K9958V90M)	

Packaging

#	Item Code	Package Description	Marketing Start Date	Marketing End Date
1	NDC:79370-001-01	5 L in 1 CONTAINER; Type 0: Not a Combination Product	10/09/2020	

Marketing Information

Marketing Category	Application Number or Monograph Citation	Marketing Start Date	Marketing End Date
OTC monograph not final	part333A	10/09/2020	

Labeler - SINBIO, S.A. DE C.V. (951576424)**Registrant** - VR INTERNATIONAL LLC (074181935)**Establishment**

Name	Address	ID/FEI	Business Operations
SINBIO, S.A. DE C.V.		951576424	pack(79370-001) , manufacture(79370-001) , label(79370-001)